

[www.ortimplant.com](http://www.ortimplant.com)

# STERILIZATION CONTAINERS

AND ACCESSORIES

" We produce for your health... "



## CONTAINER CATALOG

 **ortimplant**

[www.ortimplant.com](http://www.ortimplant.com)



"We  
produce  
for your  
health..."

Ortimplant provides a complete platform of aluminum sterilization containers suitable for all surgical sets which are convenient, hygienic, safe & which can be used without wrapping with surgical textiles and which can maintain sterilization safety for a long time

**Basic Model  
Half Size Containers**

**Basic Model  
Special Size Containers**

**Basic Model  
Middle Size Containers**

**Basic Model  
Full Size Containers**

**Bio-Barrier Model  
Half Size Containers**

**Bio-Barrier Model  
Middle Size Containers**

**Bio-Barrier Model  
Full Size Containers**

**Mini Containers**

**Dental Containers**

**Flat Containers**

# ORTIMPLANT PRODUCT GROUPS

## ● TRAUMA PRODUCTS

We produce superior designed and functional products made off ISO 5832 - 1 Stainless steel and ISO 5832 - 3 Titanium high quality raw materials, suitable with human anatomy and meticulous quality standard

## ● STERILIZATION CONTAINERS

We produce aluminum sterilization containers suitable for all surgical sets which are convenient, hygienic, safe & which can be used without wrapping with surgical textiles and which can maintain sterilization safety for a long time

## ● TRAUMA INSTRUMENT SETS

Orthopedic surgery instrument sets, including carefully designed surgical instruments aluminum boxes, silicone holders suitable for the placement of the surgical instruments for all orthopedic surgeries, provide great convenience to surgeons during surgery

## ● STERILIZATION BASKETS

We produce perforated stainless steel baskets suitable for all sets of surgical instruments that can be used in sterilization containers or alone

## ● ORTHOPEDIC SURGICAL INSTRUMENTS

We develop and produce ergonomic, convenient and functional, stainless steel orthopedic surgical instruments which increase the success of surgical procedures

## ● ORTHOPEDIC STERILIZATION CONTAINERS

Aluminum perforated boxes, screw racks, silicone holders, surgical instrument trays are produced for all surgical instrument sets, which allow the orthopedic surgical instruments to be stored, transported and sterilized by wrapping surgical textile

"With its obsessive Quality policy...."

With its obsessive quality policy and with its open structure to development, Ortimeplant goes on its constitution under today's modern institutional structure

**Endoscopy Containers**

**Implant Containers**

**Paper Filter**

**Textile Filter**

**Teflon Filter**

**Paper Filter Pack Of 100 For Mini Size Containers**

**Paper Filter Pack Of 100 For Dental Size Containers**

**Plastic Security Seal**

**Coloured Identification Label**

**Paper Label with Indicator**

**Indicator Strips**

**Silicone Mats**

Our expert team has years of experience in development of business and strategy..

From raw materials to services after sales, Ortimeplant considers to apply Total Quality Management System, With respect to this policy

**Ortimeplant has been certified its Quality with the certificates of CE, ISO 9001, ISO 13485.**



Ortimplant manufactures orthopaedics implants, sterilization containers and accessories in line with world norms with its high technologic workshop.

Ortimplant increases its success in orthopaedic field everyday by manufacturing new and developed products as a result of investigation and research studies of the specialist physicians and experienced engineers according to needs of global sterilization market.

Together with modern background and educated technical team Ortimplant aims qualified manufacturing as a goal.

From raw materials to services after sales Ortimplant considers to apply Total Quality Management. With respect to this policy Ortimplant has been certified its Quality with the certificates of CE, ISO 9001, ISO 13485.

With the aim of more rapid and continuous satisfaction of customer desires, Ortimplant educates the employees in advanced stages to increase customer satisfaction under all conditions. With periodically given education seminars Ortimplant is in the aim of reaching institutional quality from personal quality and with direction of this goal Ortimplant succeeded the communication between employees by meeting all management and technical teams together in social and cultural activities.

Ortimplant attributes to national and international sectoral fairs to actively present its current and new products to the medical sector in the world.

With its obsessive Quality policy and with its open structure to developments Ortimplant goes on its constution under today's modern instutional structure.





# INDEX

|                                                             |    |
|-------------------------------------------------------------|----|
| • Basic Model - Half Size Containers.....                   | 6  |
| • Basic Model - Special Size Containers.....                | 6  |
| • Basic Model - Middle Size Containers.....                 | 7  |
| • Basic Model - Full Size Containers.....                   | 7  |
| • Bio-Barrier Model - Half Size Containers.....             | 8  |
| • Bio-Barrier Model - Middle Size Containers.....           | 8  |
| • Bio-Barrier Model - Full Size Containers.....             | 9  |
| • Mini Containers.....                                      | 9  |
| • Endoscopy Containers.....                                 | 10 |
| • Implant Containers.....                                   | 10 |
| • Dental Containers.....                                    | 11 |
| • Flat Containers.....                                      | 11 |
| • Plastic Security Seal.....                                | 12 |
| • Coloured Identification Label.....                        | 12 |
| • Paper Label with Indicator.....                           | 12 |
| • Indicator Strips .....                                    | 12 |
| • Silicone Mats.....                                        | 12 |
| • Paper Filter .....                                        | 13 |
| • Textile Filter.....                                       | 13 |
| • Teflon Filter .....                                       | 13 |
| • Paper Filter Pack Of 100 For Mini Size Containers.....    | 13 |
| • Paper Filter Pack Of 100 For Dental Size Containers ..... | 13 |

## Basic Model - Half Size Containers

outside dimensions / mm  
**285x280x100**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O310.10               | O311.10           |
| YELLOW | O320.10               | O321.10           |
| RED    | O330.10               | O331.10           |
| BLUE   | O340.10               | O341.10           |
| GREEN  | O350.10               | O351.10           |
| BLACK  | O360.10               | O361.10           |



outside dimensions / mm  
**285x280x135**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O310.13               | O311.13           |
| YELLOW | O320.13               | O321.13           |
| RED    | O330.13               | O331.13           |
| BLUE   | O340.13               | O341.13           |
| GREEN  | O350.13               | O351.13           |
| BLACK  | O360.13               | O361.13           |

outside dimensions / mm  
**285x280x150**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O310.15               | O311.15           |
| YELLOW | O320.15               | O321.15           |
| RED    | O330.15               | O331.15           |
| BLUE   | O340.15               | O341.15           |
| GREEN  | O350.15               | O351.15           |
| BLACK  | O360.15               | O361.15           |

outside dimensions / mm  
**285x280x200**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O310.20               | O311.20           |
| YELLOW | O320.20               | O321.20           |
| RED    | O330.20               | O331.20           |
| BLUE   | O340.20               | O341.20           |
| GREEN  | O350.20               | O351.20           |
| BLACK  | O360.20               | O361.20           |


outside dimensions / mm  
**285x280x260**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O310.26               | O311.26           |
| YELLOW | O320.26               | O321.26           |
| RED    | O330.26               | O331.26           |
| BLUE   | O340.26               | O341.26           |
| GREEN  | O350.26               | O351.26           |
| BLACK  | O360.26               | O361.26           |

## Basic Model - Special Size Containers

outside dimensions / mm  
**370x280x100**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O910.10               | O911.10           |
| YELLOW | O920.10               | O921.10           |
| RED    | O930.10               | O931.10           |
| BLUE   | O940.10               | O941.10           |
| GREEN  | O950.10               | O951.10           |
| BLACK  | O960.10               | O961.10           |



outside dimensions / mm  
**370x280x135**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O910.13               | O911.13           |
| YELLOW | O920.13               | O921.13           |
| RED    | O930.13               | O931.13           |
| BLUE   | O940.13               | O941.13           |
| GREEN  | O950.13               | O951.13           |
| BLACK  | O960.13               | O961.13           |



outside dimensions / mm  
**370x280x50**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O910.15               | O911.15           |
| YELLOW | O920.15               | O921.15           |
| RED    | O930.15               | O931.15           |
| BLUE   | O940.15               | O941.15           |
| GREEN  | O950.15               | O951.15           |
| BLACK  | O960.15               | O961.15           |

## Basic Model - Middle Size Containers

outside dimensions / mm  
**465x280x100**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O210.10               | O211.10           |
| YELLOW | O220.10               | O221.10           |
| RED    | O230.10               | O231.10           |
| BLUE   | O240.10               | O241.10           |
| GREEN  | O250.10               | O251.10           |
| BLACK  | O260.10               | O261.10           |



outside dimensions / mm  
**465x280x135**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O210.13               | O211.13           |
| YELLOW | O220.13               | O221.13           |
| RED    | O230.13               | O231.13           |
| BLUE   | O240.13               | O241.13           |
| GREEN  | O250.13               | O251.13           |
| BLACK  | O260.13               | O261.13           |



outside dimensions / mm  
**465x280x150**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O210.15               | O211.15           |
| YELLOW | O220.15               | O221.15           |
| RED    | O230.15               | O231.15           |
| BLUE   | O240.15               | O241.15           |
| GREEN  | O250.15               | O251.15           |
| BLACK  | O260.15               | O261.15           |

## Basic Model - Full Size Containers

outside dimensions / mm  
**580x280x100**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O110.10               | O111.10           |
| YELLOW | O120.10               | O121.10           |
| RED    | O130.10               | O131.10           |
| BLUE   | O140.10               | O141.10           |
| GREEN  | O150.10               | O151.10           |
| BLACK  | O160.10               | O161.10           |



outside dimensions / mm  
**580x280x135**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O110.13               | O111.13           |
| YELLOW | O120.13               | O121.13           |
| RED    | O130.13               | O131.13           |
| BLUE   | O140.13               | O141.13           |
| GREEN  | O150.13               | O151.13           |
| BLACK  | O160.13               | O161.13           |

outside dimensions / mm  
**580x280x150**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O110.15               | O111.15           |
| YELLOW | O120.15               | O121.15           |
| RED    | O130.15               | O131.15           |
| BLUE   | O140.15               | O141.15           |
| GREEN  | O150.15               | O151.15           |
| BLACK  | O160.15               | O161.15           |

outside dimensions / mm  
**580x280x200**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O110.20               | O111.20           |
| YELLOW | O120.20               | O121.20           |
| RED    | O130.20               | O131.20           |
| BLUE   | O140.20               | O141.20           |
| GREEN  | O150.20               | O151.20           |
| BLACK  | O160.20               | O161.20           |

outside dimensions / mm  
**580x280x260**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O110.26               | O111.26           |
| YELLOW | O120.26               | O121.26           |
| RED    | O130.26               | O131.26           |
| BLUE   | O140.26               | O141.26           |
| GREEN  | O150.26               | O151.26           |
| BLACK  | O160.26               | O161.26           |


## Bio-Barrier Model - Half Size Containers

outside dimensions / mm  
**285x280x100**

| COLOUR | Bottom Non-Perforated |
|--------|-----------------------|
| GRAY   | O310.10B              |
| YELLOW | O320.10B              |
| RED    | O330.10B              |
| BLUE   | O340.10B              |
| GREEN  | O350.10B              |
| BLACK  | O360.10B              |



Bottom Non Perforated



outside dimensions / mm  
**285x280x150**

| COLOUR | Bottom Non-Perforated |
|--------|-----------------------|
| GRAY   | O310.15B              |
| YELLOW | O320.15B              |
| RED    | O330.15B              |
| BLUE   | O340.15B              |
| GREEN  | O350.15B              |
| BLACK  | O360.15B              |

outside dimensions / mm  
**285x280x135**

| COLOUR | Bottom Non-Perforated |
|--------|-----------------------|
| GRAY   | O310.13B              |
| YELLOW | O320.13B              |
| RED    | O330.13B              |
| BLUE   | O340.13B              |
| GREEN  | O350.13B              |
| BLACK  | O360.13B              |

outside dimensions / mm  
**285x280x200**

| COLOUR | Bottom Non-Perforated |
|--------|-----------------------|
| GRAY   | O310.20B              |
| YELLOW | O320.20B              |
| RED    | O330.20B              |
| BLUE   | O340.20B              |
| GREEN  | O350.20B              |
| BLACK  | O360.20B              |


outside dimensions / mm  
**285x280x260**

| COLOUR | Bottom Non-Perforated |
|--------|-----------------------|
| GRAY   | O310.26B              |
| YELLOW | O320.26B              |
| RED    | O330.26B              |
| BLUE   | O340.26B              |
| GREEN  | O350.26B              |
| BLACK  | O360.26B              |


## Bio-Barrier Model - Middle Size Containers

outside dimensions / mm  
**465x280x100**

| COLOUR | Bottom Non-Perforated |
|--------|-----------------------|
| GRAY   | O210.10B              |
| YELLOW | O220.10B              |
| RED    | O230.10B              |
| BLUE   | O240.10B              |
| GREEN  | O250.10B              |
| BLACK  | O260.10B              |



Bottom Non Perforated



outside dimensions / mm  
**465x280x135**

| COLOUR | Bottom Non-Perforated |
|--------|-----------------------|
| GRAY   | O210.13B              |
| YELLOW | O220.13B              |
| RED    | O230.13B              |
| BLUE   | O240.13B              |
| GREEN  | O250.13B              |
| BLACK  | O260.13B              |

outside dimensions / mm  
**465x280x150**


| COLOUR | Bottom Non-Perforated |
|--------|-----------------------|
| GRAY   | O210.15B              |
| YELLOW | O220.15B              |
| RED    | O230.15B              |
| BLUE   | O240.15B              |
| GREEN  | O250.15B              |
| BLACK  | O260.15B              |



## Bio-Barrier Model - Full Size Containers

outside dimensions / mm  
**580x280x100**

| COLOUR | Bottom Non-Perforated |
|--------|-----------------------|
| GRAY   | Ø110.10B              |
| YELLOW | Ø120.10B              |
| RED    | Ø130.10B              |
| BLUE   | Ø140.10B              |
| GREEN  | Ø150.10B              |
| BLACK  | Ø160.10B              |



outside dimensions / mm  
**580x280x135**

| COLOUR | Bottom Non-Perforated |
|--------|-----------------------|
| GRAY   | Ø110.13B              |
| YELLOW | Ø120.13B              |
| RED    | Ø130.13B              |
| BLUE   | Ø140.13B              |
| GREEN  | Ø150.13B              |
| BLACK  | Ø160.13B              |

Bottom Non Perforated



outside dimensions / mm  
**580x280x150**

| COLOUR | Bottom Non-Perforated |
|--------|-----------------------|
| GRAY   | Ø110.15B              |
| YELLOW | Ø120.15B              |
| RED    | Ø130.15B              |
| BLUE   | Ø140.15B              |
| GREEN  | Ø150.15B              |
| BLACK  | Ø160.15B              |

outside dimensions / mm  
**580x280x200**

| COLOUR | Bottom Non-Perforated |
|--------|-----------------------|
| GRAY   | Ø110.20B              |
| YELLOW | Ø120.20B              |
| RED    | Ø130.20B              |
| BLUE   | Ø140.20B              |
| GREEN  | Ø150.20B              |
| BLACK  | Ø160.20B              |


outside dimensions / mm  
**580x280x260**

| COLOUR | Bottom Non-Perforated |
|--------|-----------------------|
| GRAY   | Ø110.26B              |
| YELLOW | Ø120.26B              |
| RED    | Ø130.26B              |
| BLUE   | Ø140.26B              |
| GREEN  | Ø150.26B              |
| BLACK  | Ø160.26B              |


## Mini Containers

outside dimensions / mm  
**300x140x50**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | Ø410.04               | Ø411.04           |
| YELLOW | Ø420.04               | Ø421.04           |
| RED    | Ø430.04               | Ø431.04           |
| BLUE   | Ø440.04               | Ø441.04           |
| GREEN  | Ø450.04               | Ø451.04           |
| BLACK  | Ø460.04               | Ø461.04           |



Bottom Non Perforated



Bottom Perforated



outside dimensions / mm  
**300x140x70**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | Ø410.07               | Ø411.07           |
| YELLOW | Ø420.07               | Ø421.07           |
| RED    | Ø430.07               | Ø431.07           |
| BLUE   | Ø440.07               | Ø441.07           |
| GREEN  | Ø450.07               | Ø451.07           |
| BLACK  | Ø460.07               | Ø461.07           |

outside dimensions / mm  
**300x140x100**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | Ø410.10               | Ø411.10           |
| YELLOW | Ø420.10               | Ø421.10           |
| RED    | Ø430.10               | Ø431.10           |
| BLUE   | Ø440.10               | Ø441.10           |
| GREEN  | Ø450.10               | Ø451.10           |
| BLACK  | Ø460.10               | Ø461.10           |

Endoscopy Containers

outside dimensions / mm

**250x60x50**

| COLOUR | Bottom Perforated |
|--------|-------------------|
| GRAY   | O711.05           |
| YELLOW | O721.05           |
| RED    | O731.05           |
| BLUE   | O741.05           |
| GREEN  | O751.05           |
| BLACK  | O761.05           |



Bottom  
Perforated ( Non Filter )



outside dimensions / mm

**450x70x70**


| COLOUR | Bottom Perforated |
|--------|-------------------|
| GRAY   | O711.07           |
| YELLOW | O721.07           |
| RED    | O731.07           |
| BLUE   | O741.07           |
| GREEN  | O751.07           |
| BLACK  | O761.07           |

Implant Containers

outside dimensions / mm

**310x138x65**

| COLOUR | Bottom Perforated |
|--------|-------------------|
| GRAY   | O611.06           |
| YELLOW | O621.06           |
| RED    | O631.06           |
| BLUE   | O641.06           |
| GREEN  | O651.06           |
| BLACK  | O661.06           |



Bottom  
Perforated ( Non Filter )



outside dimensions / mm

**500x155x75**

| COLOUR | Bottom Perforated |
|--------|-------------------|
| GRAY   | O611.07           |
| YELLOW | O621.07           |
| RED    | O631.07           |
| BLUE   | O641.07           |
| GREEN  | O651.07           |
| BLACK  | O661.07           |

## Dental Containers

outside dimensions / mm  
**310x190x50**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O510.04               | O511.04           |
| YELLOW | O520.04               | O521.04           |
| RED    | O530.04               | O531.04           |
| BLUE   | O540.04               | O541.04           |
| GREEN  | O550.04               | O551.04           |
| BLACK  | O560.04               | O561.04           |



outside dimensions / mm  
**310x190x65**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O510.06               | O511.06           |
| YELLOW | O520.06               | O521.06           |
| RED    | O530.06               | O531.06           |
| BLUE   | O540.06               | O541.06           |
| GREEN  | O550.06               | O551.06           |
| BLACK  | O560.06               | O561.06           |

outside dimensions / mm  
**310x190x100**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O510.10               | O511.10           |
| YELLOW | O520.10               | O521.10           |
| RED    | O530.10               | O531.10           |
| BLUE   | O540.10               | O541.10           |
| GREEN  | O550.10               | O551.10           |
| BLACK  | O560.10               | O561.10           |

Bottom Non Perforated



Bottom Perforated




outside dimensions / mm  
**310x190x130**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O510.13               | O511.13           |
| YELLOW | O520.13               | O521.13           |
| RED    | O530.13               | O531.13           |
| BLUE   | O540.13               | O541.13           |
| GREEN  | O550.13               | O551.13           |
| BLACK  | O560.13               | O561.13           |


## Flat Containers

outside dimensions / mm  
**285x280x55**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O810.05               | O811.05           |
| YELLOW | O820.05               | O821.05           |
| RED    | O830.05               | O831.05           |
| BLUE   | O840.05               | O841.05           |
| GREEN  | O850.05               | O851.05           |
| BLACK  | O860.05               | O861.05           |



Bottom Non Perforated



Bottom Perforated



outside dimensions / mm  
**285x280x95**

| COLOUR | Bottom Non-Perforated | Bottom Perforated |
|--------|-----------------------|-------------------|
| GRAY   | O810.08               | O811.08           |
| YELLOW | O820.08               | O821.08           |
| RED    | O830.08               | O831.08           |
| BLUE   | O840.08               | O841.08           |
| GREEN  | O850.08               | O851.08           |
| BLACK  | O860.08               | O861.08           |

## Plastic Security Seal


| CAT.NO. | PRODUCT NAME                      |
|---------|-----------------------------------|
| O011.01 | Pack Of 100 Plastic Security Seal |



## Coloured Identification Label


### Coloured Identification Labes

| COLOUR | Bottom Non-Perforated |
|--------|-----------------------|
| GRAY   | O020.01               |
| YELLOW | O020.02               |
| RED    | O020.03               |
| BLUE   | O020.04               |
| GREEN  | O020.05               |
| BLACK  | O020.06               |




## Paper Label With Indicator

| CAT.NO. | PRODUCT NAME                                                              |
|---------|---------------------------------------------------------------------------|
| O010.01 | Pack Of 100 Paper Label With Indicator For Half-Midi-Full Size Containers |
| O010.02 | Pack Of 100 Paper Label With Indicator For Dental and Mini Containers     |



## Indicator Strips

| CAT.NO. | PRODUCT NAME                         |
|---------|--------------------------------------|
| O813.00 | Pack Of 250 Class 6 Indicator Strips |




## Paper Filter Pack Of 100 For Mini Size Containers



| CAT.NO. | PRODUCT NAME                                                 |
|---------|--------------------------------------------------------------|
| O012.02 | Paper Filter Pack Of 100 For Mini Size and Dental Containers |


## Silicone Mats



| MM      | CAT.NO. |
|---------|---------|
| 180x280 | O050.01 |
| 245x245 | O050.03 |
| 270x125 | O050.04 |

| MM      | CAT.NO. |
|---------|---------|
| 380x230 | O050.05 |
| 440x230 | O050.06 |
| 520x230 | O050.07 |

## Teflon Filter



| CAT.NO. | PRODUCT NAME                                                   |
|---------|----------------------------------------------------------------|
| O014.01 | PTFE Filter For Half, Special, Middle and Full Size Containers |
| O014.02 | PTFE Filter For Mini Containers                                |
| O014.03 | PTFE Filter For Dental Containers                              |

## Textile Filter



| CAT.NO. | PRODUCT NAME                                                             |
|---------|--------------------------------------------------------------------------|
| O013.01 | Single Textile Filter For Half, Special, Middle and Full Size Containers |
| O013.02 | Single Textile Filter For Mini Containers                                |
| O013.03 | Single Textile Filter For Dental Containers                              |


## Paper Filter




| CAT.NO. | PRODUCT NAME                                                                |
|---------|-----------------------------------------------------------------------------|
| O012.01 | Paper Filter Pack Of 100 For Half, Special, Middle and Full Size Containers |



# PRODUCTGROUPS



TRAUMA PRODUCTS




TRAUMA INSTRUMENT SETS



ORTHOPEDIC SURGICAL INSTRUMENTS



STERILIZATION CONTAINERS



STERILIZATION BASKETS




ORTHOPEDIC STERILIZATION CONTAINERS





**Ortimplant**  
**manufactures orthopaedics implants,**  
sterilization containers and accessories in line with world norms  
with its high technologic workshop.



# ortimplant

[www.ortimplant.com](http://www.ortimplant.com)

Ortimplant Samsun  
Kirazlık Mahallesi, 1034. Cadde No:14 Tekkek y  
SAMSUN / TURKEY  
[export@ortimplant.com](mailto:export@ortimplant.com)

Workshop  
Organize Sanayi Sitesi Vali Mehmet ErdoĖan Cebeci  
Bulvarı No:34 Tekkek y / SAMSUN / TURKEY

Tel : +90 362 266 97 28  
Fax : +90 362 266 58 72  
[info@ortimplant.com](mailto:info@ortimplant.com)

“  
**Quality  
in  
Health**